



Recombinant *Saccharomyces cerevisiae* Topoisomerase I damage affected protein 7 (TDA7)

Product Code	CSB-CF518968SVR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	E7QJS8
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Zymaflore VL3) (Baker's yeast)
Sequence	MNSNSTIGRTTLGESDTISLSFSEPSSSLNSRSTDVVFASTSTLVPQQGSLTSL PPVSST ATPTYYSTSLTYDETLHTSIDVSSTSTLVSSSTDSSSSSEQDTYSSQYDPATSSY SIITPS MSIFSSTSPMSSSSSITSEWSSLTSTPTLSSSATSLSSSWSSLSSPSSLLVSS SLSLSL SSSYSDTKLFSFDRSSIFSPSTPTVISPSYTYLSSISATSFQISTTSELSSSWFS TISS PSTTSNKDTTFPSSSRNTSTSFYSSSLSTNDFSTISKSSKLSPSASSSTVSI SIST SVPT SSSVSSSSSKVPSNRPSSSSSSDDTTSAYSSTYTFQSLQSTTSSSIPPTTQTPS TSTIST SPIPTSSQVFNTVAISSSEDSKIYYFYTQTYDITDSSTTFVTGLPTTIAVAKSEV TSFS APSSITADMSFYQHWLDGSLDNNKNQGTSKTNTGTIVGSSVGGILICVLV VWFMLV RKRKAKRHFKEKENDSFCHEIGRRTGFPTTAQAKEASLQAQDSGSQQRNTETAS ANNPFSNE FNFKARGNPPPVPVPPRNVAMNGSFQNMRSNFMQENRFSYGSSFTYSSLG SSTQGGFST LSSNSIRLGRGLDNDISHDERNTVQNN SQGFLREII
Source	in vitro E.coli expression system
Target Names	TDA7
Protein Names	Recommended name: Topoisomerase I damage affected protein 7
Expression Region	1-636
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.