



# Recombinant Danio rerio LETM1 and EF-hand domain-containing protein 1, mitochondrial (letm1)

<b>Product Code</b>	CSB-CF638067DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q1LY46
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Sequence</b>	VHTTGRLDDSKVERSLRTLKDRNKKLEEGGPVYSPTVDAEPVRRRTIRQRVID EVKHYYH GFRLLWIDTTIAVRMLWRVLNGHILSRERRRQFLRTCADVFRLLPFLVFIIVPFM EFLLP VALKLFPNMLPSTFETQSKKEERLKKELRVKLEMAKFLQDTIEEIALRNKASKG NVTEEF STFFQKIRDSGEIPSNEQIIRFSKLFEDLTLDNLTRPQLVALCKLLELQSIGTNN FLRF QLIMKLRAIRADDKLI AEEGVDSL TANELQAACRVRGMRALGVTEERLREQLK QWLELHL NQHIPTSLLLLSRAMFLPDTLSPADQLKTTLQNLPEIMAKEAQVKVAELDFSKV DNKTKL ETTLQEEAAIRQENRERELERLADAAEKAKEQTQSQEAEVLEVEGERRVDAEH ALSSVDV AIHSETLRDTAPVLEGIKGEEITKEEIDMLSDACTKLKEQKNLLTKEKEELEDLK DDVQE YSEDL EEIKRELSKTGQEKAVEESKASQRLSKRVNRMIGRMDKIITELEKDKMV LDGQMD SETTPPIGENLISINELITVMKQIQNIPEHKLLSIADALDENKDGKIDDDVIKVVELI D KEDIDISTNQVAEIMVMLQKEEKLMEKEKAKEKVEKEQAQKIQN
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	letm1
<b>Protein Names</b>	Recommended name: LETM1 and EF-hand domain-containing protein 1, mitochondrial Alternative name(s): Leucine zipper-EF-hand-containing transmembrane protein 1
<b>Expression Region</b>	114-757
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.