



Recombinant Mouse Sulphydryl oxidase 2 (Qsox2)

Product Code	CSB-CF659699MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3TMX7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	VGPREGGGARLYREGSDAVWLLDSGVSRSATGNSSAAWLQFHSSWCGHCI GYAPTWRAL AADVRDWAAAIRVAALDCAEEKNQDVCRTYDIHFYPTFRYFKAFTKEFTTGEN FKGPDRE LRTVRQTMIDFLQNHTEGTWPPACPLDPIQSSDILSFMDSHSGQYHAIVFESN GSYVGR EVILDLIPYENIMVSRALDTDKAFLGTLGITSVPSCYLIYPNGSHGLVNVAKPLRS FFSS HLKSLPDVRKKSLFLPEKSNKEEKSEVVVWKEFDRAKLYTADLESGLHYLLRV ELAAHRS LAGAQLKTRDFVTVVAKLFPGRPAVKKLLLETLQEWLANLPLDKIPYNAILDLVN NKMQUI SGIFLTSHVKWVGCQGSRLRLRGYPCSLWKLFTLTQASTHPEALAGTGFE GHPQAVLQ AIRRYIRTFFGCKEKGHEFEEMAKESMDSVKTPDQAVLWLWRKHNMVNSRLA GHLSEDPK FPKVPWPPTDLCACHEEIKGLDSWNEGQVLLFLKQHYSRDNLVDAYSVDQG SPGEWEAQ GREQEEGKGLNPSGKSWRHHDGTSLRPPHILGPRTDLSKSLHHRLDLRLQSP QGPQALKE AKAVVPFLGVGFSSLDMSLCVVLYVASSFLMIMYFFFRVRSKRWKVRLYHPA V
Source	in vitro E.coli expression system
Target Names	Qsox2
Protein Names	Recommended name: Sulphydryl oxidase 2 EC= 1.8.3.2 Alternative name(s): Quiescin Q6-like protein 1
Expression Region	39-692
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.