



Recombinant Human Protein SERAC1 (SERAC1)

Product Code	CSB-CF853443HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96JX3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MSLAAYCVICRRIGTSTSPPKSGTHWRDIRNIIKFTGSLILGGSLFLTYEVLALK KAVT LDTQVVEREKMKSYYVHTVSLDKGENHGIAWQARKELHKAVRKVLATSakilR NPFADP FSTVDIEDHECAVWLLLRKSKSDDKTTRLEAVREMSETHHWHDYQYRIIAQAC DPKTLIG LARSEESDLRFFLLPPPLPSLKEDSSTEEELRQLLASLPQTElDECiQYFTSLAL SESSQ SLAAQKGGLWCFGGNGLPYAESFGEVPSATVEMFCLEAIVKHSEISTHCDKIE ANGLLQL LQRLYRLHKDCPKVQRNIMRVIGNMALNEHLHSSIVRSGWVSIMAEAMKSPHI MESSHAA RILANLDRETVQEKYQDGVYVLHPQYRTSQQPIKADVLFIHGLMGAAFKTWRQQ DSEQAVI EKPMEDEDRYTTCWPKTWLAKDCPALRIISVEYDTSLSDWRARCPMERKSIAF RSNELLR KLRAAGVGDRPVVWISHSMGGLLVKKMLLEASTKPEMSTVINNTRGIIFYSVPH HGSRLA EYSVNIRYLLFPSLEVKELSKDSPALKTLQDDFLEFAKDKNFQVLNFVETLPTYI GSMIK LHVVPVESADLGIGDLIPVDVNHNLNICKPKKKDAFLYQRTLQFIREALAKDLEN</p>
Source	in vitro E.coli expression system
Target Names	SERAC1
Protein Names	Recommended name: Protein SERAC1
Expression Region	1-654
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.