



Recombinant Danio rerio Protein SERAC1 (serac1)

Product Code	CSB-CF716322DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5SNQ7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSVSALRLIRVRRLSTSGPAVKRALPWRDIRRIAKVTGAIVFGGCVFITYEVVTL NQALT IDTSAILQEKQKSYIYPHSTNREQESLASGLTIKTRRELHKAARKFLEITSRVLH HPLD EHLSHLDADPHECALWVLLKRSRSADRAVRHLAVQELAHNHHWRDYQYQTA AQVVDQRTA VALARIPNVDLRFFLPPPPLPHTEDDISIEDGLRQLLASLPQSDVDQCVQYFTSL ALRES SQSLASQRGGLWCFGGNGLPYAQSLTSTPSEKVETFCQLQALVQHSKVRSHCE HIVANGGL QLLQRVYQLRRDSPKIQRNIVRIIGNLALNENLHTTIVQSGWMSVLAEMIQSPHI MQASH AARALANLDRDAVRQKYQDGVYILHPQCRTNQPIKADVLFVHGLLGAAFKTR QKDCDVT DDEKLEGVREDYTECWPKSWLAADCPNLRILSVEYDTHLSDWNSKCPVENQR KSLAFRSQ ELLRKLKDAGVGERPVIWVAHSMGGLLVKKMLLDAAKDPDLSSLIKNTKGILFY SVPHHG TFMAEYSVSVRYLLFPSIEVKELCRDSPALRDLNENFLNIAKDREFKVLSFAETV PTYIG PMLKILVPAHSADLGIGDLIQVDVDHLNICKPEKKDFTLYKRTLQFIQDALGGR RIK
Source	in vitro E.coli expression system
Target Names	serac1
Protein Names	Recommended name: Protein SERAC1
Expression Region	1-658
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.