



Recombinant Human Syntabulin (SYBU)

Product Code	CSB-CF865166HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NX95
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGPLRESKKEHRVQHHDKEISRSRIPRLILRPHMPQQQHKVSPASESPFSEEE SREFNPS SSGRSARTVSSNSFCSDDTGCPSSQSVSPVKTPSDAGNSPIGFCPGSDEGFT RKKCTIGM VGEGSIQSSRYKKESKSGLVKPGSEADFSSSSSTGSISAPEVHMSTAGSKRSS SSRNRGP HGRSNGASSHKPGSSPSPREKDLLSMLCRNQLSPVNIHPSYAPSSPSSNS GSYKGSDC SPIMRRSGRYMSCGENHGVPPNPEQYLTPQQKEVTVRHLKTKLKESERRL HERESEIV ELKSQ LARMREDWIEEECHRVEAQLALKEARKEIKQLKQVIETMRSSLADKDK GIQKYFV DINIQNKKLESLLQSMEMAHSGSLRDELCLDFPCDSPEKSLTLNPPLDTMADG LSLEEQV TGEGADRELLVGDSIANSTDLFDEIVTATTTESGDLELVHSTPGANVLELLPIVM GQEEG SVVVERAVQTDVVPYSPAISELIQSVLQKLQDPCPSSLASPDESEPD SMESFP ELSALV VDLTPRNPNSAILLSPVETPYANVDAEVHANRLMRELDFAACVEERLDGVIPLA RGGVVR QYWSSSFLVDLLAVAAPVVPTVLWAFSTQRRGGTDPVYNIGALLRGCCVVALHS LRRTAFR IKT
Source	in vitro E.coli expression system
Target Names	SYBU
Protein Names	Recommended name: Syntabulin Alternative name(s): Golgi-localized syntaphilin-related protein Syntaxin-1-binding protein
Expression Region	1-663
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.