



Recombinant Mouse Syntabulin (Sybu)

Product Code	CSB-CF805899MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BHS8
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MGPLRESKKEQRVQHQEKEISRSRIPRLILRPHRPQQQQQQQNKVSPASESP FSEESRE FNPSSSGRSARTISSNSFCSDDTGCPSSQSVSPVKTPSDTGHSPIGFCPGSDE DFTRKKC RIGMVGEGSIQSARHKKEPKGGIIPGSEADFSSSSSTGSISAPEVHMSTTGNK RASFSR NRGPHGRSNGASSHKSGSSPPSPREKDLVSM LCRNPLSPSNIHPSYAPSSPS SSNSGSYK GSDCSPVMRRSGRYMSCGENHGVKPPNPEQYLTP LQQKEVTVRHLR TKLKE SERLHERE SEIMELKSQLARMREDWIEEECHRVEAQLALKEARKEIKQLKQVIETMRSSLAD KDKGIQ KYFVDINIQNKKLESLLQSMEMAHNSSLRDELCLDFSFDSPEKSLPLSSTFDKL PDGLSL EEQITEEGADSELLVGDSMAEGTDLLDEMVTATTTESSGLEFVHSTPGPQALK ALPLVSH EEGIAVMEQAVQTDVVPFSPAISELIQSVLKLQDYCPTSSASPDESGADSMESF SESISA LMLDLTPRSPNSAILLSPVEIPFSKGAMEAHANRLMRELDFAAYTEERLDSVLS LSQGSV VRQYWSSNFLVDLLAVAAPVVPTVLWAFSTQRGGTDPVYNIGALLRGCCVVA LHSLRRTA FHMKT</p>
Source	in vitro E.coli expression system
Target Names	Sybu
Protein Names	Recommended name: Syntabulin Alternative name(s): Golgi-localized syntaphilin-related protein Short name= m-Golsyn Syntaxin-1-binding protein
Expression Region	1-665
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.