



# Recombinant Arabidopsis thaliana Probable beta-1,3-galactosyltransferase 18 (B3GALT18)

<b>Product Code</b>	CSB-CF840428DOA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q8RX55
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Sequence</b>	<p>MKKPKLSKVEKIDKIDLFSSLWKQRSVRVIMAIGFLYLVIVSVEIPLVFKSWSSSS  VPLD  ALSRLEKLNNEQEPQVEIIPNPPLEPVSYVSNPTIVTRTDLVQNKVREHHRGV  LSSLRF  DSETFDPSSKDGSELHKSAAKEAWQLGRKLWKELESGRLEKLVEKPEKNKPD  SCPHSVSL  TGSEFMNRENKLMELPCGLTLGSHITLVGRPRKAHPKEGDWSKLVSQFVIELQ  GLKTVEG  EDPPRILHFNPRLKGDWSKPKVIEQNSCYRMQWGPAQRCEGWKSRDDEETV  DSHKCEKW  IRDDNYSEGSRARWWLNRLIGRRKRVKVEWPFVVEEKLFVLTLSAGLEGY  HINVDGKH  VTSFPYRTGFTLEDATGLTVNGDIDVHSVFVSLPTSHPSFAPQRHLELSKRW  QAPVVPD  GPVEIFIGILSAGNHFSERMAVRKSWMQHVLITSAKVVARFFVALHGRKEVNVE  LKKEAE  YFGDIVLPYMDSYDLVVLKTVAICEHGALAFSAKYIMKCDDDTFVKLGAVINEV  KKVPE  GRSLYIGNMNYHKLPLRGGKWAVTYEEWPEEDYPPYANGPGYVLSSDIARFI  VDKFERHK  LRLFKMEDVSVGMMVEHFKNTTNPVDYRHSLRFCQFGCVENYYTAHYQSPR  QMICLWDKL LRQNKPECCNMR</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	GALT5
<b>Protein Names</b>	Recommended name: Probable beta-1,3-galactosyltransferase 18 EC= 2.4.1.-
<b>Expression Region</b>	1-672
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.