



Recombinant Human Sulfhydryl oxidase 2 (QSOX2)

Product Code	CSB-CF757950HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6ZRP7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	ARRSPPPRAARLPRLLVLLAAAVGPGAGGAARLYRAGEDAVWVLDSGSVRG ATANSSAA WLVQFYSSWCGHCIGYAPTWRALAGDVRDWASAIRVAALDCMEEKNQAVCH DYDIHFYPT FRYFKAFTKEFTTGENFKGPDRELRTVRQTMIDFLQNHTEGSRPPACPRLDPI QPSDVLS LLDNRGSHYVAIVFESNSSYLGREVILDLIPYESIVVTRALDGDKAFLEKLGVS VPSCY LIYPNGSHGLINVVKPLRAFFSSYLKSLPDVRKKSLPLPEKPHKEENSEIVVWR EFDKSK LYTVDLESGLHYLLRVELAAHKS LAGAELKTLKDFVTVLAKLFPGRPPVKKLLE MLQEWL ASLPLDRIPYNAVLDLVNNKMRISGIFLTNHIKWVGCQGSRSELRGYP CSLWKL FHTLTV EASTHPDALVGTGFEDDPQAVLQTMRRYVHTFFGCKEKGEGHFEEMAKESMD SVKTPDQAI LWLWKKHNMVNGRLAGHLS EDPFRFPKLQWPTDLC PACHEEIKGLASWDEG HVLTF LKQH YGRDNLLDTYSADQGDSSEG GTLARGE EEEEEKRLTPPEVSHGDRDTQSVRPP GALGPRPAL PESLHHS LDGKLQSLDGP GAHKEVGGAAPFLGVDFSSLDMSLCVVLYVASSLF LMVMYFF FRVRSRRWKVKHHHPAV
Source	in vitro E.coli expression system
Target Names	QSOX2
Protein Names	Recommended name: Sulfhydryl oxidase 2 EC= 1.8.3.2 Alternative name(s): Neuroblastoma-derived sulfhydryl oxidase Quiescin Q6-like protein 1
Expression Region	22-698
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.