



# Recombinant *Arabidopsis thaliana* Probable beta-1,3-galactosyltransferase 20 (B3GALT20)

<b>Product Code</b>	CSB-CF416687DOA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A7XDQ9
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Arabidopsis thaliana</i> (Mouse-ear cress)
<b>Sequence</b>	<p>MKRVKSESFRGVYSSRRFKLSHFLLAIAGFYLVFLAFKFPHFIEMVAMLSGDTG LDGALS DTSLDVSLSGSLRNDMLNRKLEDEDHQSGPSTTQKVSPEEKINGSKQIQPLLF RYGRISG EVMRRRNRTIHMSPFERMADEAWILGSKAWEDVDKFEVDKINESASIFEGKVE SCPSQIS MNGDDLKANKRIMLLPCGLAAGSSITILGTPQYAHKESVPQRSRLTRSVMVL VSQFMVE LQGLKTGDGEYPPKILHLNPRKGDWNHRPVIEHNTCYRMQWGV AQRC DGTP SKKDADVL VDGFRRCWKWTQNDIIDMVDSKESKTTSWFKRFIGREQKPEVTWSFPFAEGK VFVLT LRA GIDGFHINVGGRHVSSFPYRPGFTIEDATGLAVTGDVDIHSIHATSLSTSHPSFS PQKAI EFSSEWKAPPLPGTPFRLFMGVLSATNHFSERMAVRKTWMQHPSIKSSDVVA RFFVALNP RKEVNAMLKKEAEYFGDIVILPFMDRYELVVLKTIACEFGVQNV TAPYIMKCDD DTFIR VESILKQIDGV SPEKSLYMGNLNLRHRPLRTGKWTVTWEEWPEAVYPPYANG PGYIISSN IAKYIVSQNSRHKLR LFKMEDVSMGLWVEQFNASMQPVEYSHSWKFCQYGCT LNYYTAHY QSPSQMMCLWDNLLKGRPQCCNFR</p>
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	GALT2
<b>Protein Names</b>	Recommended name: Probable beta-1,3-galactosyltransferase 20 EC= 2.4.1.-
<b>Expression Region</b>	1-684
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.