



Recombinant Rhizobium radiobacter Protein virD4 (virD4)

Product Code	CSB-CF362802RKV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09817
Product Type	Transmembrane Protein
Immunogen Species	Rhizobium radiobacter (Agrobacterium tumefaciens) (Agrobacterium radiobacter)
Sequence	<p>MNSSKTSPQRMTLSIVCSLAAGFCAASCYVTFRRGFNGEAIMTFDVFAFWYET PLYLGYA STVFWRGLSVVIFTSLAVLSSQLIISLRNQKHHGTARWAEIGEMQHAGYLQRYS RIKGPI FGKTCGPLWFGSYLTNGEQPHSLVVAPTRAGKGVGIVIPTLLTFKGSVIALDVK GELFEL TSRARKASGDAVFKFSPLDSEKTHCYNPVLDIATLPPERQFTETRRLAANLIT AKGKGA EGFIDGARDLFVAGILTCIERGTPTIGAVYDLFAQPGEKYKLFQAQAEESLNKEA QRIFD NMAGNDTKILTSYTSVLGDGGLNLWADPLIKAATSRSDFSVYDLRRKKTICIYLC VSPNDL EVLAPLMRLLFQQLVSILQRSLPGKDECHEVLFLLDEFKHLGKLEAIETAITTIAG YKGR FMFIIQSLSALTGTYYDDAGKQNFLSNTGVQVFMATADDETPTYISKAIGEYTFKA RSTSY SQASMFHDHNIQISDQGAALLRPEQVRLDDQSEIVLIKGRPPLKLRKVQYYSDR TLKRLF ERQMGSLPEPAPLMLSDYSNDQVQSHLAEIANFNEDAAPRNRTVAEDHGSVK VGADIPER VMGINGDEDQADAREIPPELGCASRIDAGSGRSTAIVGPDYCTSAKIEVRTGTA CEMKLN WVQFVSGGMLRFDKVYARSQLSLYEVS</p>
Source	in vitro E.coli expression system
Target Names	virD4
Protein Names	Recommended name: Protein virD4
Expression Region	1-688
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.