



Recombinant Danio rerio UPF0668 protein C10orf76 homolog (zgc:63733)

Product Code	CSB-CF747484DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6PGW3
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MAQTEKKAGLLRRTSSSKKPLKEKVVLMYDEIFMKEDPTKNSPRFWDELFLMK VNLEYLE GKLESLDGDEFMKIKDNINSLFHHCVQALAEEHQIKVVNALQTLCALFRGVHQK NKSATG FDIINMLIGFDKAELRMKELMESLDSLLCGDGSESLKSLCLKLLLCLVTVDNIS QNTIL EYVMINSIFEAILQILSDVSSRSQHGYDSVLLALLVNYRKYESVNPYIVKLSIVD DEPT LDGLGMVIHHALTEYNRQYKEKEEENQGGFFSTLSSMVGSMFIADADEKLSVQ TNEAILL ALYEAVHLNRNFITVLAQSHPEIDMAATPATPVPASPTTPLGTTTPSLDLMNPN ELPLDP NLQTSNLLITFLKYSSIVMQDTKDEHRLNSARLCLILTCIAEDQYADAFHDDN MFRV NLHRMPMRHRKKAADKNIPCRPLVCAVLDLMVEFVVTHMMKEFPMDLYVRCV QIIHKLIC YQKKCRVRLHYTWRELWSALINLLKFLLSNETGLLAKHNIFQLALQVVNLFNMF ITYGDT FLPTPSSYDELYYEIIRMHVFDNLYCMVLRVSTNAGQWKEPASKVTHALVNV RAIINH NPKIESYA AVNHISQLSEDQVLEVRSNYDTLTLKLQDGLDQFERYSQPREA GFFKELVRSISLNRKNVSLTTMSQDVLLKEFSSIS</p>
Source	in vitro E.coli expression system
Target Names	zgc:63733
Protein Names	Recommended name: UPF0668 protein C10orf76 homolog
Expression Region	1-689
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.