



Recombinant *Schizosaccharomyces pombe* Uncharacterized transcriptional regulatory protein PB1A11.04c (SPAPB1A11.04c)

| | |
|--------------------------|---|
| Product Code | CSB-CF867210SXV |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q9HDX1 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) |
| Sequence | <p>MIAAVEKDAVPRKRRVISQVSQACIRCRQKKIKCSGEEKPSCQACSNNKVECIW PDRPNMR GKKPSSAARTLSPLTGIVQPYLGSHPFSPTILSPSSPNMDFTYINSFGDYQKQ LRLPTS RKLLQLWDIFLKTSYTELFGIFNKAQITSQIASDTAPPVLVLCICAHAAARFSQEEV NRFK SSTAASDYANQAFSLLPCRFDISLTNITCLLLLCLIELGSCRGAKAWLLLGMA LRMVD SLDLGNEINDNPLTMGSNTVSWTEAEQKRRVYWACFFVERLLSTGYMAPSKL RSLSLSLQ KTSIQLPCPEPNFLFNQPILTELDGSLPENTQSDTTSMAPYQRSLQLEFMTGS LIRLSN LWSEISRWALCGGYTKDITPPWLNQSQFHRSFELLKAWHENLPPRAVWSYTN YSAYSSPG ESAGACYTFMHLLYHTTLTYLLRNVLDFPEKSRQKSKLFSSVSQRFGQPPT VWMDMIL DQVITSADFITKLSKDPLNYIMSPFVGFSLTAATIHMLLKFCVVNIDQNYISSRL VHV DHQILQDRSKYWKINQAMLVTLQRLYNFYRFQYLEEQSLYNFKIPGFPLCILEY GIVEDQ SMKSNPDLNSFSKELFSRGTNELLNNGDDTQSGNSTPAMGVSEIRTDILDEE VPVDPLI TSILDDGRWWEEMFGSERKAGFKETVFEDMNGRSIRL</p> |
| Source | in vitro E.coli expression system |
| Target Names | SPAPB1A11.04c |
| Protein Names | Recommended name: Uncharacterized transcriptional regulatory protein PB1A11.04c |
| Expression Region | 1-697 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.