



# Recombinant Human DnaJ homolog subfamily C member 14 (DNAJC14)

<b>Product Code</b>	CSB-CF744334HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q6Y2X3
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	MAQKHPGERGLYGAHHSGGASLRTLGPSVDPEIPSFSGLRDSAGTAPNGTRC LTEHSGPK HTQHNPNAHWLDP SHGPPGGPGPPRDAEDPDQSETSSEEEESGVDQELSKEN ETGNQK DGN SFLSIPSACNCQGTPGIPEGYPYSEGGNGSSSNFCHHCTSPALGEDELEEEYDD EESLKFP SDFSRVSSGKKPPSRQRHRFPTKEDTREGGRRDPRSPGRHRLGRKRSQAD KRKGLGLWG AEELCQLGQAGFWWLIELLVLVGEYVETCGHLIYACRQLKSSDLDLFRVWMG VWTGRLGG WAQVMFQFLSQGFYCGVGLFTRFLKLLGALLLALALFLGFLQLGWFLVGLG DRLGWRD KATWLF SWLDSPALQRCLTLRDSRPWQRLVRIVQGWLELPWVKQINIRQ GNAPVASGR YCQPEEEVARLLTMAGVPEDELNPFHVLGVEATASDVELKKAYRQLAVMVHP DKNHHPRA EEAFKVLRAAWDIVSNAEKRKEYEMKRMAENELSRVNEFLSKLQDDLKEAM NTMMCSRC QGKHRRFEMDREPKSARYCAECNRLHPAEEGDFWAESSMLGLKITYFALMD GKVYDITEW AGCQRVGISPDTHRVPYHISFGSRIPGTRGRQRATPDAPPADLQDFLSRIFQV PPGQMPN GNFFAAPQPAPGAAAASKPNSTVPKGAEKPKRRKKVRRPFQR
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	DNAJC14
<b>Protein Names</b>	Recommended name: DnaJ homolog subfamily C member 14 Alternative name(s): DnaJ protein homolog 3 Dopamine receptor-interacting protein of 78 kDa Short name= DRIP78 Human DnaJ protein 3 Short name= hDj-3
<b>Expression Region</b>	1-702
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.