



Recombinant Escherichia coli O139:H28 Bifunctional protein aas (aas)

Product Code	CSB-CF426104EJD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZQU2
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Sequence	<p>MLFSFFRNLCRVLYRVRVTGDTQALKGERVLITPNHVSFIDGILLGLFLPVRPVF AVYTS ISQQWYMRWLKSFIDFVPLDPTQPMAIKHLVRLVEQGRPVVIFPEGRITTTGSL MKIYDG AGFVAAKSGATVIPVRIEGAELTHFSRLKGLVKRRLFPQITLHILPPTQVAMPDA PRARD RRKIAGEMLHQIMMEARMAVRPRETLYESLLSAMYRFGAGKKCVEDVNFTPD SYRKLLTK TLFVGRILEKYSVEGERIGLMLPNAGISAAVIFGAIARRRIPAMMNYTAGVKGLT SAITA AEIKTIFTSRQFLDKGKLWHLPEQLTQVRWVYLEDLKADVTTADKVVIFAHLLM PRLAQV KQQPEEEALILFTSGSEGHKPGVVHSHKSILANVEQIKTIADFTTNDRFMSALPL FHSFG LTVGLFTPLLTGAEVFLYPSPLHYRIVPELVYDRSCTVLFGTSTFLGHYARFANP YDFYR LRYVVAGAEKLQESTKQLWQDKFGLRILEGYGVTECAPVVSINVPMAAKPGTV GRILPGM DARLLSVPGIEEGGRLQLKGPNIMNGYLRVEKPGVLEVPTAENVRGEMERGW YDTGDIVR FDEQGFVQIQGRAKRFKAIAGEMVSLEMVEQLALGVSPDKVHATAIKSDASKG EALVLFT TDNELTRDKLQQYAREHGVPELAVPRDIRYLKQMPLLGSGKPDFVTLKSWVD EAEQHDE</p>
Source	in vitro E.coli expression system
Target Names	aas
Protein Names	<p>Recommended name: Bifunctional protein aas Including the following 2 domains: 2-acylglycerophosphoethanolamine acyltransferase EC= 2.3.1.40 Alternative name(s): 2-acyl-GPE acyltransferase Acyl-[acyl-carrier-protein]- phospholipi</p>
Expression Region	1-719
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.