



Recombinant Tyrosine-protein kinase wzc (wzc)

Product Code	CSB-CF834796EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8X7L9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MTEKVKQHAAPVTGSDEIDIGRLVGTVIEARWWVIGITAVFALCAVVYFFFATPI YSADA LVQIEQSSGNSLVQDIGSALANKPPASDAEIQILRSRLVLGKTVDDLDLDIAVSK NTFPI FGAGWDRMLMGRQNETVKVTTFNRPKEMEDQVFTLNVLNKNYTLSSDGGFS ARGQAGQIL KKEGVTLMVEAIHARPGSEFTVTKYSTLGMINQLQNSLTVTENGKDAGVLSLTY TGEDRE QIRDILNSIARNYQEQNIERKSAEASKSLAFLAQQQLPEVRSRLDVAENKLNFR QDKDSV DLPLEAKAVLDSMVNIDAQLNELTFKEAEISKLYTKVHPAYRTLLEKRQALEDEK AKLNG RVTAMPKTQQEIVRLTRDVESGQQVYMQLLNKEQELKITEASTVGDVRIVDPAI TQPGVL KPKKGLIILGAILGLMLSIVGVLLRSLFNRGIESPQVLEEHGISVYASIPLSEWQK ARD SVKTIKGIKRYKQSQLLAVGNPTDLAIEAIRSLRSLHFAMMQAQNVLMMTG VSPSIGK TFVCANLAAVISQTNKRVLIDCDMRKGYTHELLGTNNVNLSEILIGQGDITTA AKPTS IAKFDLIPRGQVPPNPSELLMSERFAELVNWASKNYDLVLIDTPPILAVTDAAIV GRHVG TTLMVARYAVNTLKEVETSLSRFEQNGIPVKGVILNSIFRRASAYQDYGYEYE YKSDAK
Source	in vitro E.coli expression system
Target Names	wzc
Protein Names	Recommended name: Tyrosine-protein kinase wzc EC= 2.7.10.-
Expression Region	1-720
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.