



# Recombinant Escherichia coli Tyrosine-protein kinase wzc (wzc)

<b>Product Code</b>	CSB-CF303499ENV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P76387
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Sequence</b>	<p>MTEKVKQHAAPVTGSDEIDIGRLVGTVIEARWWVIGITTVFALCAVVYFFFATPI            YSADA            LVQIEQNSGNSLVQDIGSALANKPPASDAEIQLRSLVLGKTVDDLDDIAVSK            NTFPI            FGAGWDRMLMGRQNETVKVTTFNRPKEMADQVFTLNVLDNKNYTLSSDGGFS            ARGQAGQML            KKEGVTLMVEAIHASPGEFTVTKYSTLGMINQLQNSLTVTENGKDAGVLSLTY            TGEDRE            QIRDILNSIARNYQEQNIERKSAEASKSLAFLAQLPEVRSRLDVAENKLNFR            QDKDSV            DLPLEAKAVLDSMVNIDAQLNELTFKEAEISKLYTKVHPAYRTLLEKRQALEDEK            AKLNG            RVTAMPKTQQEIVRLTRDVESGQQVYMQLLNKEQELKITEASTVGDVRIVDPAI            TQPGVL            KPKKGLIILGAILGLMLSIVGVLLRSLFNRGIESPQVLEEKGISVYASIPLSEWQK            ARD            SVKTIKGIKRYKQSLLAVGNPTDLAIEAIRSLRTSLHFAMMQAQNVLMMTG            VSPSIGK            TFVCANLAAVISQTNKRVLIDCDMRKGYTHELLGTNNVNGLSEILIGQGDITTA            AKPTS            IAKFDLIPRGQVPPNPSELLMSERFAELVNWASKNYDLVLIDTPPILAVTDAIV            GRHVG            TTLMVARYAVNTLKEVETSLSRFEQNGIPVKGVILNSIFRRASAYQDYGYEYE            YKSDAK</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	wzc
<b>Protein Names</b>	Recommended name: Tyrosine-protein kinase wzc EC= 2.7.10.-
<b>Expression Region</b>	1-720
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.