



Recombinant *Saccharomyces cerevisiae* Mitochondrial outer membrane protein YKR018C (YKR018C)

Product Code	CSB-CF336130SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36114
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	<p>MFKVFGFGAKEEIPELSQEEKTKAILKQAHDFEQALRAMDYVLDDNADEGLAL LDESDAK EASDQTINALARGVIEFLEATLGFEAEEMKKASATLAKAEALSLKSRERAQKIGL KSSSL YPPGTVYAVTYTESCLLHALLMIFSESMMEAAKAILKLRKSYMLQEILETIKAA NKAKK LKITSGSEDKESTPATFITGGDAFNSVDIPYELTPEEQKDKDLLQFAEQIHSMRT ERLSG AHIGNSPAINRLRGELGLQAMEDLPEEEITDHKVLSDDIDLSQATIDEFVHSGVN LCFGI LQVVISLLPPAIGAVLSVVGFRGSREEGLRLVWKATKQRNVHGCIGLLALMFYY DGPFQF TDDDFDIPAAVKDSSNSEDESEEMDGPTLLHPGKILEDALLQSRALFPNSAL WLLNEAR MLSGKGRLEEAVALMDSIDVSKIRMRQVKSLMIFDRAITLIHLHQYDRAAEDILS LLDIS DWSHAFYTYFAGCCYLENWRMCEMGLMKSDKKDEYQKKAELIFTSVNLLGK KTFKSKNL PLDRFILRKVEQFKAKKEELGVENPLDGIATSPVHEIAYFYNGYNRMSEEHLEL TKKMLT EYRNPAIEALDSDQELIKDLLVSLTLRRLGHIQEGCDILDEKVLPKFFSIQNGKVK YIKK TEDPWAYPTALYERALFTWKLEGMDGLPESKEWLLRAQGYADDYELSTRVG MKIKAAIDR VDHSL</p>
Source	in vitro E.coli expression system
Target Names	YKR018C
Protein Names	Recommended name: Mitochondrial outer membrane protein YKR018C
Expression Region	1-725
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.

**Protein Length**

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.