



# Recombinant Mouse SUN domain-containing protein 2 (Sun2)

<b>Product Code</b>	CSB-CF804816MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q8BJS4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	<p>MSRRSQRLTRYSQDDNDGGSSSSGASSVAGSQGTVFKDSPLRTLKRKSSNM  KHLSPAPQL  GPSSDSHTSYYESVRESYIGSPRAVSLARSALLDDHLHSEPYWSGDLRGR  RRRGTGGS  ESSKANGLTAESKASEDFFGSSSGYSSDDLGYTDSQHQSSGSRRLRSAASR  AGSFVWTL  VTFPGRLFGLLYWWIGTTWYRLTTAASLLDVFVLTLSRHFSLNLKSFLWFLLLL  LLLTGL  TYGAWHFYPLGLQTLQPAVVSWWAAKESRKQPEVWESRDASQHFQAEQRV  LSRVHSLERR  LEALADFSSNWQKEAIRLERLELRQGAAGHGGGSSLSHEDALSLEGLVSRR  EATLKED  LRRDTVAHIQEELATLRAEHHQDSEDLFKKIVQASQESEARVQQLKTEWKSMT  QEAFAQES  SVKELGRLEAQLASLRQELAALTLKQNSVADEVGLLPQKIQAARADVESQFPD  WIRQFLL  GDRGARSGLLQRDEMHAQLQELENKILTKMAEMQGKSAREAAASLGQILQKE  GIVGVTEE  QVHRIVKQALQRYSEDRIGMVDYALESGGASVISTRCSSETYETKTALLSLFGIPL  WYHSQ  SPRVILQPDVHPGNCWAFQGPQGFAVVRLSARIRPTAVTLEHVPKALSPNSTI  SSAPKDF  AIFGFDEDLQQEGTLLGTFAYDQDGEPIQTFYFQASKMATYQVVELRILTNGW  HPEYTCI YRFRVHGEPAH</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Sun2
<b>Protein Names</b>	Recommended name: SUN domain-containing protein 2 Alternative name(s): Protein unc-84 homolog B Sad1/unc-84 protein-like 2
<b>Expression Region</b>	1-731
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.