



Recombinant Arabidopsis thaliana Receptor protein kinase-like protein At4g34220 (At4g34220)

Product Code	CSB-CF850092DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q94C77
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	LNTDGVLLLTFKYSILTDPLSVLRNWNYYDDATPCLWTGVCTELGKPNTPDMF RVTSLVL PNKHLGSLITPDLFSIPYLRILDSSNFFNGSLPDSVFNATELQSSISLGSNNLSGD LPKS VNSVTNLQLLNLSANAFTGEIPLNISLLKNLTVVSLSKNTFSGDIPSGFEAAQILD LSSN LLNGSLPKDLGGKSLHYLNLSHNKVLGEISPNFAEKFPANATVDLSFNLTGPI PSSLSL LNQKAESFSGNQELCGKPLKILCSIPSTLSNPPNISETTSPAIAVKPRSTAPINPL TEKP NQTGKSKLKPSTIAAITVADIVGLAFIGLLVLYVYQVRKRRRYPESKFSFFKFC LEKNE AKKSKPSTTEVTVPEspeAKTTCGSCIILTGGRYDETSTSESdVENQQTvQAFT RTDGGQ LKQSSQTQLVTVdGETRLDLDTLLKASAYILGTTGTGIVYKAVLENGTAFavRRI ETESC AAAKPKEFEREVRAIAKLRHPNLVRIRGFCWGDDEKLLISDYVPNGSLLCFFTA TKASSS SSSSSSLQNPLTFEARLKIARGMARGLSYINEKKQVHGNIKPNNILLNAENEPIIT DLGL DRLMTPARESHTTGPTSSSPYQPPEWSTSLKPNPKWDVYSFGVILLELLTSKV FSVDHDI DQFSNLSDSAAEENGRFLRLIDGAIRSDVARHEDAAMACFRLGIECVSSLPQK RPSMKEL VQVLEKICVLV
Source	in vitro E.coli expression system
Target Names	At4g34220
Protein Names	Recommended name: Receptor protein kinase-like protein At4g34220
Expression Region	27-757
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.