



# Recombinant *Drosophila melanogaster* Putative mitochondrial inner membrane protein (CG6455)

<b>Product Code</b>	CSB-CF308600DLU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P91928
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Sequence</b>	<p>MYRLAVRDQCKCALQRTLQQTANNRQFGGSSSSGSGGREQQRRQQEEQQG            QGDQGYQGYQ            SLPPHMREAGFGKVVLVVSPLAAVGGVITYAKYDDDFRKLVEKNVPGAGSVIK            VALQEEP            PFKGITKNVNDQIDKVKSGIETVTSTVDSVTSKVTGLFGGGSGDDKSKKSKVEP            VKATPA            EEKRPSKPSEVSKTEAKPVSKPAAAAAPAPAAKPKDNPLPRDVVELEKAIELSA            QLAVKE            YNVAIGVLKGFNDDVRKVVDKAVENGENSELWTTLNKRASARDTAVATAERAA            REAQEKIV            ACEIALSAAATAQNAKKVEAVRDKIKKLVDHIGNVKDELYRHKDTASVSDKYW            RNVEKAR            NYFIDEIESIFPGLSLADKKNLSKEDLDLFIHAYTHVLAYQKELQRLQTDGELR            LKRA            IDSVRGDNDSEALRAQLEYHLEAERRKLAVERNQKKIFHIHAESDKLLRLQLKKQ            AEAHAD            HIKDIVAQRETDLTRSFKRELEDKATEKANYKLQLAGMLGKLRGMDAALAER            ADAERTA            NQAQALWAACQALWASVRAATPGVHYKDRLRPLKNEINAIKVAKGDDLVA            VLESVPKE            AQERGVYPEDALRERFLNVERVARRLALVPEEGAGLPIYFLSYLQSLFILRPDN            PISKDE            LENKPFDYSKLDTYDILNRARYHVDRSDFLQALKYMNLLQGASRKIAGEWMKE            ARLMLET QQAANTLMAHAAASGLLYL</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Mitofilin
<b>Protein Names</b>	Recommended name: Putative mitochondrial inner membrane protein
<b>Expression Region</b>	1-739
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.