



Recombinant Mouse Translocation protein SEC63 homolog (Sec63)

Product Code	CSB-CF837736MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8VHE0
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	AGQQFQYDDSGNTFFYFLTSFVGLIVIPATYYLWPRDQNAEQIRLKNIRKVYGR CMWYRL RLLKPQPNIIPTVKKIVLLAGWALFLFLAYKVSKTDREYQEYNPYEVLNLDPGAT VAEIK KQYRLLSLKYHPDKGGDEVFMRIAKAYAALTDEESRKNWEEFGNPDGPQAT SFGIALPA WIVDQKNSILVLLVYGLAFMVILPVVVGSWWYRSIRYSGDQILIRTTQIYTYFVY KTRNM DMKRLIMVLGASEFDPQYNKDSTSRPTDNILIPQLIREIGSINLKKNEPPLTCPY SLKA RVLLLSHLARMKIPETLEEDQQFMLKKCPALLQEMVNVICQLIIMARSREEREF RAPTLA SLENCMKLSQMAVQGLQQFKSPLLQLPHIEEDNLRRVSNHKKYKIKTIQDLVSL KESDRH SLLHFLEDEKYEEVMAVLGSFPYVTMDIKSQVLDDSDSNNITVGS�VTVLVKLT RQTMAE VFEKEQSICAAEEQPTEDGQSDANKIKAKGGWQQKNKGPKKMPKSKKKKPLK KKPTTVPL PQAKQQKQKQANGVVGSEAAIKEEEDDISDKGSDSEEEETNRDSQSEKEDGS DRESDREQ DEKQSKDDEAEWQELQQSIQRKERALLETKSKITHPVYSLYFPEEKQEWWWL YIADRKEQ TLISMPYHVCTLKDTEEVELKFPAPGKPGNYQYTVFLRSDSYMGLDQIKPLKLE VHEAKP VPENHPQWDTAIEGDEDQEDSEGFEDSFEEEEEEEEEGGD
Source	in vitro E.coli expression system
Target Names	Sec63
Protein Names	Recommended name: Translocation protein SEC63 homolog
Expression Region	2-760
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.