



# Recombinant Pongo abelii Translocation protein SEC63 homolog (SEC63)

<b>Product Code</b>	CSB-CF725287PYX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q5R660
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Sequence</b>	AGQQFQYDDSGNTFFYFLTSFVGLIVIPATYYLWPRDQNAEQIRLKNIRKVYGR CMWYRL RLLKPQPNIIPTVKKIVLLAGWALFLFLAYKVSKTREYQEYNPYEVLNLDPGAT VAEIK KQYRLLSLKYHPDKGGDEVMFMRIAKAYAALTDEESRKNWEEFGNPDGPQAT SFGIALPA WIVDQKNSILVLLVYGLAFMVILPVVVGSWWYRSIRYSGDQILIRTTQIYTYFVY KTRNM DMKRLIMVLGASEFDPQYNKDATS RPTDNILIPQLIREIGSINLKKNEPPLTCPY SLKA RVLLLSHLARMKIPETLEEDQQFMLKKCPALLQEMVNVICQLIVMARNREEREF RAPTLA SLENCMKLSQMAVQGLQQFKSPLLQLPHIEEDNLRRVSNHKKYKIKTIQDLVSL KESDRH TLLHFLEDEKYEEVMAVLGSPYVTMDIKSQVLDDSDSNNITVGS�VTVLVKLT RQTMAE VFEKEQSICAAEEQPAEDGQGETNKNRTKGGWQQKSKGPKKTAKS KKKKKP LKKKPTPVL LPQSKQQKQKQANGVVGNEAAVKEDEEEVSDKGS DSEEEETNRDSQSEKDD GSDRDSRE QDEKQNKDDEAEWQELQQSIQRKERALLETKSKITHPVYSLYFPEEKQEWWW LYIADRKE QTLISMPYHVCTLKDTEEVELKFPAPGKPGNYQYTVFLRSDSYMGLDQIKPLKL EVHEAK PVPENHPQWDTAIEGDEDQEDSEGFEDSFEEEEEEEEEDDD
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SEC63
<b>Protein Names</b>	Recommended name: Translocation protein SEC63 homolog
<b>Expression Region</b>	2-761
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.