



Recombinant Transmembrane GTPase fzo-1 (fzo-1)

Product Code	CSB-CF637567CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q23424
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	<p>MSGTASLVHTLPASGDSNHRGLHSLKNSRRAADNEPLLRFREAKKVLGDVYG ELKDNVAE LEGVYKDIKENDFVSSEQREEIEAIGDSIKTIMDTFQRDNMKVVFFGRTSNGKS TTINAM LHEKVLPQGMGHTTCCFLQVEGSEGEVGHLLQDDNPQKIDMKMLGKIGHALS DENSDLPA MGQDSLLKVFHPKSESSEGCERLLQNDVVILDSPGVDLSPDFDSWIDKHCLDAD VFVLVSN AESTLTQAEKNFFLRVAKKLSKPNVFILNNRWDASAAETENIEDVKKQHLTRFR QFLVDE LEVCSEREVNDRIFFVSSREVLESRLKARGLVQKAYQAEGHGTRALEFQNFER HFEHCIS RSAIHTKFEAHNRRRAHEMIGKMRLNLSVLTSAAEQRSKLQNNLNESTRTFNE CRVNFTQ FEKAYREQTEQLRAEVHLKVSADFFEEIARLDAIIDRFEQPFDGSSSGMTKYKE DLAIFV DKCLSSDLEARCTGGLMSRIWNLENDMFQYVTKILAEPYQNKLEEVWRYRAP FKFSICVD VPALVNDFHEDLEFRFTFGLHAIIRRIIAYRSGQPVTAINLLTPLSLKQQSEKN SVRD AEASAASEEQAMMTQMVLTSAAFLANGSLGVLVGGIVYKAVGWRVIAVGGA AYAGLYAW ERMWNSGAKEQHLKEQFRSHLAARMQQVSTAHTHHCETQAIREMDQVFDG LKATVGGVH REMKNDLDVQKTQIDAVDSTIRTLGTIKGKAVFLLRNLEQFASSYLRS DSPPTP</p>
Source	in vitro E.coli expression system
Target Names	fzo-1
Protein Names	Recommended name: Transmembrane GTPase fzo-1 EC= 3.6.5.-
Expression Region	1-774
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.