



Recombinant Human Metal transporter CNNM4 (CNNM4)

Product Code	CSB-CF764727HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6P4Q7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPVGGGGRPVGGPARGRLLLAAPVLLVLLWALGARGQGSPQQGTIVGMRL ASCNKSCGT NPDGIIFVSEGSTVNLRLYGYSLGNISSNLISFTEVDDAETLHKSTSCLELTKDLV VQQL VNVSRGNTSGVLVLTkFLRRSESMKLYALCTRAQPDGPWLKWTDKDSLLFM VEEPGRFL PLWLHILLITVLLVLSGIFSGNLGLMALDPMELRIVQNCGTEKERRYARKIEPIR RKGN YLLCSLLLGNVLVNTSLTILLDNLIGSGLMAVASSTIGIVIFGEILPQALCSRHGLA VGA NTILLTKFFMLLTFPLSFPISKLLDFFLGQEIRTVYNREKLMEMLKVTEPYNDLVK EELN MIQGALELRKTVEDIMTQLQDCFMIRSDAILDFNTMSEIMESGYTRIPVFEDEQ SNIVD ILYVKDLAFVDPDDCTPLKTITRFYNHPVHFVFHDTKLDAMLEEFKKGKSHLAIV QKVNN EGEGDPFYEVGLVLTLEDVIEEIIKSEILDESDMYTDNRSRKRVSSEKNKRDFSAF KDADN ELKVKISPQLLLAAHRFLATEVSQFSPSLISEKILLRLLKYPDVIQELKFDEHNKY YARH YLYTRNKPADYFILILQGKVEVEAGKENMKFETGAFSYYGTMALTSVPSDRSP AHPTPLS RSASLSYPDRTDVSTAATLAGSSNQFGSSVLGQYISDFSVRALVDLQYIKITRQ QYQNGL LASRMENSPQFPIDGCTTHMENLAEKSELPVDETTLNERNSSLHKASHEN AI
Source	in vitro E.coli expression system
Target Names	CNNM4
Protein Names	Recommended name: Metal transporter CNNM4 Alternative name(s): Ancient conserved domain-containing protein 4 Cyclin-M4
Expression Region	1-775
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.