



# Recombinant Pan troglodytes Reticulon-1 (RTN1)

<b>Product Code</b>	CSB-CF689275EQV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q5IS59
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Pan troglodytes (Chimpanzee)
<b>Sequence</b>	<p>MAAPGDPQDELLPLAGPGSQWLRDRGEGEDEAVTPKGATPAPQAGEPSPGL  GARAREAAS  REAGSGPARQSPVAMETASTGVAGVSSAMDHTFSTTSKDGECSYTSLSIDI  CYPPQEDS  TYFTGILQRENGHVTISESPEELGTPGSSLPDVPGIESRGLFSSDSGIEMTPAE  STEVNK  ILADPLDQMKAAYKYIDITRPEEVKHQEQNHPELEDKDLDFKNKDTDISIKPEG  VREPD  EPAPVEGKIIKDHLLEESTFAPYIDDLSEEQRRAPQITTPVKITLTEIEPSVETTTQ  EKT  PEKQDICKPSPDTVPTVTVSEPEDDSPGSITPPSSGTEPSAAESQGKGSISED  ELITAI  KEAKGLSYETAESPRPVGQLADRPEVKARSGPPTIPSPLDHEASSAESGDSEI  ELVSEDP  MAAEDALPSGYVSFGHVGPPSPASPSIQYSILREEREAELDSELIIESCDAS  SASEES  PKREQDSPPMKPGALDAIREETGVRAEERAPSRRLAEPASFLDYPSTEPQP  GPELPPGD  GALEPETPTLPRKPEEDASSHQSPAATKGGPGLGPGAPPPLLFLNKQKAIDLLY  WRDIKQ  TGIVFGSFLLLLFSLTQFSVSVVAYLALAALSATISFRIYKSVLQAVQKTDEGHP  FKAY  LELEITLSQEQIQKYTDCLQFYVNSTLRELRLFLVQDLVDSLKFAVLMWLLTYV  GALFN  GLTLLLMAVVSMFTLPVVVYVKHQAQIDQYLGLVRTHINAVVAKIQAKIPGAKRH  AE</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	RTN1
<b>Protein Names</b>	Recommended name: Reticulon-1
<b>Expression Region</b>	1-776
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.