



Recombinant Human Putative protein cTAGE-6 (CTAGE6P)

Product Code	CSB-CF771435HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q86UF2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEEPGATPQPYLGLVLEELRRVVAALPESMRPDENPYGFPSELVVCAAVIGFF VVLLFLW RSFRSVRSRLYVGREQKLGATLSGLIEEKCKLLEKFSLIQKEYEGYEVESLED ASFEKA AAEEARSLEATCEKLNRSNSELEDEILCLEKDLKEEKSQDELMADISKSI QSLEDE SKSLKSQIAEAKIICKTFKMSEERRAIAIKDALNENSQLQTSHKQLFQQEAEVWK GEVSE LNKQKITFEDSKVHAEQVLNDKENHIKTLTGHLPMMDQAAVLEEDTTDDDNL ELEVNSQ WENGANLDDPLKGALKKLIHAAKLNVSLSLEGERNHIIIQLSEVDKTKEELTEH IKNLQ TQQESLQSENIYFESENQKLQQKLKIMTEFYQEDEMPLYRKLTVREENYRIEEEE KLSKVE EKLSRATEQLETYRKLAKDLEEEELERTVHFYQKQVISYEKRGHDNWLAARTAE RNLSDLR KENAHNKQKLTETELKFELLEKDPNALDVSNTAFGREHAPNGPAPLGQRSSET RAFLSPQ TLLEDPLGLSPVLPEGGGRGPRGPGNPLDHQITNERGEPSCDRLTDPHRAPS DTGSLSSP VEQDCKMMFPPPGQSYPDALPPQREDRFYSNSERLSGSAEPRSFKMTSLD KMDGSMPSE MESSRNDKDDLGNLNVPDSSLPAENEATGPGFIPPLAPVRGPLFPVDTRGP FMRRGPP FPPPPPGTMFGASRGYFPPRDFPGPPHAPFAMRNIYPPRGLPPYFHPRPGFY PNPAF
Source	in vitro E.coli expression system
Target Names	CTAGE6
Protein Names	Recommended name: Putative protein cTAGE-6 Alternative name(s): cTAGE family member 6
Expression Region	1-777
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.