



Recombinant *Schizosaccharomyces pombe* Uncharacterized membrane protein C6F6.13c (SPAC6F6.13c)

Product Code	CSB-CF515660SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14244
Product Type	Transmembrane Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Sequence	<p> MPISSPGTRCSSLKDPDLQYSAESVSTEQSLGTFEESKGSITENYVQDSSV DEHDDGN WQPMEVISLEPTHLINDIDDDNEIIEEKKETEKVEESELEPRYTRVFRDEDDDQ KHQLDS EAIKLLDIADHGNEEISMDSQLEITGNILSETEKMAYAGVCRLILKLVDMVKIACFT TLPW YRGECKAALEDTIMWADKTTSCIYEHLGVTVEEQMIENLHKHSVQIDDLKIL VSAHRA QTVSSLDVAVLVDEVESDSSLSSLGKEKPVQIDVRWTVLCDLFLVLISKSLYDCR SRSLLM AVGEVLDINEFDVAKFEKHIVETIQIDDTGELEAGSSANTEAVMKLRRKVSRRK KYILMG LAGIGGGLVIGLSSGLLAPIISAGIGAAFTTVGLSGVATSGFLAGGGSAALITAGG AISG AHIGTTGMAHRKADVKTFFRPLHAQRRANVIVTVSGWMLSKEDDVRLSFATL DPIVGD YSVFWEPEMLASAGQTMNILATEVVTQSLQQVLGSTVLVSLMGALQWPLILTK LGYLIDN PWNNSLDRAKATGQLLADMLCYRSLGVRPVTLVGYSLGARVIYYCLRELEKKK EFSIEN VYLFGTPVIFKRTSWLKAASVVSGRFVNGYKKNWILGYLFRATSGGIGRVAG LRQIDCI PGIENIDVTNLVSGHLAYRESMPILLAAVGFVLDVEEVDLVSEPIPEPLRERQSQ LLYEI EAEECQNKQKELIEKSLMQGRSLSPKKSNAFFDSKKIREELKKVKKKYGSSSF NSRWY </p>
Source	in vitro E.coli expression system
Target Names	SPAC6F6.13c
Protein Names	Recommended name: Uncharacterized membrane protein C6F6.13c
Expression Region	1-778
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.



Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>