



Recombinant Human Vezatin (VEZT)

Product Code	CSB-CF872540HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9HBM0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MTPEFDEEVVFENSPLYQYLQDLGHTDFEICSSLSPKTEKCTTEGQQKPPTRV LPKQGIL LKVAETIKSWIFFSQC�KKD DLLHKLDIGFRLDSLHTILQQEVLLQEDVELIELLD PSIL SAGQSQQQENGHLPTLCSLATPNIWDL SMLFAFISLLVMLPTWWIVSSWLWV GVILFVYL VIRALRLWRTAKLQVTLKKYSVHLEDMATNSRAFTNLVRKALRLIQETEVISRG FTLVSA ACPFNKAGQHPSQHILIGLRKAVYRTL RANFQAARLATLYMLKNYPLNSESDNV TNYICVV PFKELGLGLSEEQISEEEAHNFTDGFSLPALKVLFQLWVAQSSEFFRRLALLS TANSP GPLLTPALLPHRILSDVTQGLPHAHSACLEELKRSYEFYRYFETQHQSVPQCLS KTQQKS RELNNVHTAVRSLQLHLKALLNEVILEDELEKLVCTKETQELVSEAYPILEQKLK LIQP HVQASNNCWEEAISQVDKLLRRNTDKKGKPEIACENPHCTVPLKQPTLHIAD KDPIPEE QELEAYVDDIDIDSDFRKDDFYYSQEDKERQKREHEESKRVLQELKSVLGFK ASEAERQ KWKQLLFSDHAVLKSLSPVDPVEPISNSEPSMNSDMGKVSKN DTEESNKSA TTDNEISR TEYLCENSLEGKNKDSSNEVFPQGAEERM CYQCESEDEPQADGSGLT TAP PTPRDSLQP SIKQRLARLQLSPDFTFTAGLAAEVAARSLSFTTMQEQTFGGEEEEQII EENKN EIEEK</p>
Source	in vitro E.coli expression system
Target Names	VEZT
Protein Names	Recommended name: Vezatin
Expression Region	1-779
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.