



Recombinant *Schizosaccharomyces pombe* Uncharacterized zinc metalloprotease C354.09c (SPBC354.09c)

Product Code	CSB-CF526117SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43023
Product Type	Transmembrane Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Sequence	<p>MTDEKHYVPPPKDPPSYEEVALHSALNNSAPPNDGEQNETSMEEEMIEPPS EDSSRFP LLRTKLAAIHEGWESACHSFEIRFASTFHRIPFQFLYLAVIATVILASYGYFDG VPAW RSVHHYGEDVLLNYIKGCDISDTRQQVMTLSSIPHLAGTVGDSSLLQMIMNRLY YEKGTI VDFREFYAYLNFQQLVSLSIDGDDSFHPSLIESYQVGGYDGVSIPTPATFGGSP SGFVNA PLVYANRGRIEDFEWLNSGIYVESSIVLVRANQSDFALATANAKEYNASAILIF EDTYL TSLDNLNQVYPAGPYPSANSLYRGSVANHYYYVGDPLTPGWSAHEETNRISP KDANVLP IVSIPITFNDGIELLKRLQGHGHLVKDSNWCQDLAPVLSEVWTGSKISSPGLEV NVLQDI EDKQKIINIMAQIDGYESDQILVVGAPRDSWCTGASDSSVGTSLIDVISTFANM AQDLS WKPRRTIVFASWDARQFNAIGSTEFLEYWKESLEAKAVAYINVDVAVSGDTFT ARTVPGL KKVIQRAFDVANEEMKAANIITDDFDYTSDLTSFLTIFAGIPVVNLA FERNEEN PTPMP FLGSCEDTVSWIDTFGSEYWENAARLGKIWSYLILFLANDPVVPPYDLEDEINGV GEMLKR IPEIPGANALDLRKINEEFSELLESLIRFEDEIREWKSLMMHNSYTVSVKKHPEL EGYNA KLARFERSFLDEAGLPGHEWYKHLIYGPNLRNSHSQLFPSIFDALLYGDVEAA QKEVKRI ALALDRAHNEIRFA</p>
Source	in vitro E.coli expression system
Target Names	SPBC354.09c
Protein Names	Recommended name: Uncharacterized zinc metalloprotease C354.09c EC=3.4.-.-
Expression Region	1-794
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at



4°C for up to one week.

Tag Info

Tag type will be determined during the manufacturing process.

Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.