



Recombinant Human Leucine-rich repeat-containing protein 8E (LRRC8E)

Product Code	CSB-CF743982HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6NSJ5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MIPVAEFKQFTEQQPAFKVLKPWWVDVLAEYLTVAMLMIGVFGCTLQVTQDKIIC LPNHEL QENLSEAPCQQLLPRGIPEQIGALQEYKGLKNNLDLQQYSFINQLCYETALHW YAKYFPY LVVIHTLIFMVCTSFWFKFPGTSSKIEHFISILGKCFDSPWTRALSEVSGENQK GPAAT ERAAATIVAMAGTGPGKAGEGEKEKVLAEPEKVVTEPPVVTLLDKKEGEQAKA LFEKVKK FRMHVEEGDILYTMIRQTVLKVCKFLAILVYNLVYVEKISFLVACRVETSEVTG YASFC CNHTKAHLFSKLAFCYISFVCIYGLTCIYTLWLFHRPLKEYSFRSVREETGMG DIPDVK NDFAFMLHLIDQYDSLYSKRFAVFLSEVSESRLKQLNLNHEWTPEKLRQKLQR NAAGRLE LALCMLPGLPDTVFELSEVESLRLEAICDITFPPGLSQLVHLQELSLHSPARLP FSLQV FLRDHLKVMRVKCEELREVPLWVFGLRGLEELHLEGLFPQELARAATLESLRE LKQLKVL SLRSNAGKVPASVTDVAGHLQRLSLHNDGARLVALNSLKKLAALRELELVACG LERIPHA VFSLGALQELDLKDNHLRSIEEILSFQHCRKLVTLRLWHNQAIVPEHVRKLR LEQLYL SYNKLETLPSQLGLCSGLRLLDVSHNGLHSLPPEVGLLQNLQHLALSYNALAL PEELFF CRKLRTLLLGDNQLSQLSPHVGALRALSRLELKGNRLEALPEELGNCGGLKKA GLLEDVDT LYQGLPAEVRDKMEEE</p>
Source	in vitro E.coli expression system
Target Names	LRRC8E
Protein Names	Recommended name: Leucine-rich repeat-containing protein 8E
Expression Region	1-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.

**Protein Length**

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.