



Recombinant Quaternary ammonium compound-resistance protein sugE (sugE)

Product Code	CSB-CF527874DTO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O69279
Product Type	Transmembrane Protein
Immunogen Species	Citrobacter freundii
Sequence	MSWIVLLIAGLLEVVWAIGLKYTHGFTRLTPSIITIAAMIVSIAMLSWAMRTLPGV TAYA VWTGIGAVGAAITGILLGEGASPARLLSLGLIVAGIIGLKLSTH
Source	in vitro E.coli expression system
Target Names	sugE
Protein Names	Recommended name: Quaternary ammonium compound-resistance protein sugE
Expression Region	1-105
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.