



Recombinant Rhizobium etli UPF0060 membrane protein RHE_CH01408 (RHE_CH01408)

Product Code	CSB-CF645583RKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2KAC5
Product Type	Transmembrane Protein
Immunogen Species	Rhizobium etli (strain CFN 42 / ATCC 51251)
Sequence	MTYIIYAFAAVFEIAGCFWAWLKLEKPAWWLAPGMISLALFAWLLTLVPSDA AGRTFA AYGGIYILASLSWLWLIAGRVPDRYDIGGGLICLAGASVILFAPRA
Source	in vitro E.coli expression system
Target Names	RHE_CH01408
Protein Names	Recommended name: UPF0060 membrane protein RHE_CH01408
Expression Region	1-106
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.