



Recombinant Escherichia coli O45:K1 UPF0060 membrane protein ynfA (ynfA)

Product Code	CSB-CF486880EOK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7M9T5
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Sequence	MIKTTLLFFATALCEIIGCFLPWLWLKRNASIWLLLPAIGSLALFVWLLTLHPAAS GRVY AAYGGVYVCTALIWLRLVVDGKLTLYDWTGALIALCGMLIIVAGWGRT
Source	in vitro E.coli expression system
Target Names	ynfA
Protein Names	Recommended name: UPF0060 membrane protein ynfA
Expression Region	1-108
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.