



Recombinant UPF0060 membrane protein YnfA (ynfA)

Product Code	CSB-CF812561EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8FHC6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MIKTTLLFFATALCEIIGCFLPWLWLKRNASIWLLLPA GISLALFVWLLTLHPAAS GRVY AAYGGVYVCTALIWL RVVDG VKLTLYDWTGALIALCGMLIIVAGWGRT
Source	in vitro E.coli expression system
Target Names	ynfA
Protein Names	Recommended name: UPF0060 membrane protein YnfA
Expression Region	1-108
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.