



Recombinant Enterobacter sp. Spermidine export protein MdtI (mdtI)

Product Code	CSB-CF395341EIZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4WA58
Product Type	Transmembrane Protein
Immunogen Species	Enterobacter sp. (strain 638)
Sequence	MQQFEWIHAAWLAFAIVLEIIANVFLKFSDGFRRKWFGLLSIAAVLGAFSALSQA VKGID LSVAYALWGGFGIAATLAAGWVLFQRLNRKGWIGLVLLLAGMVMIKLA
Source	in vitro E.coli expression system
Target Names	mdtI
Protein Names	Recommended name: Spermidine export protein MdtI
Expression Region	1-109
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.