



Recombinant Escherichia coli O8 Spermidine export protein MdtI (mdtI)

Product Code	CSB-CF483365E00
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LZZ1
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O8 (strain IA11)
Sequence	MAQFEWVHAAWLALAIIVLEIVANVFLKFSDGFRRKIFGLLSLAAVLAAFSALSQ AVKGID LSVAYALWGGFGIAATLAAGWILFGQRLNRKGWIGLVLLLAGMIMVKLA
Source	in vitro E.coli expression system
Target Names	mdtI
Protein Names	Recommended name: Spermidine export protein MdtI
Expression Region	1-109
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.