



Recombinant Quaternary ammonium compound-resistance protein qacE (qacE)

Product Code	CSB-CF366192ENL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGC9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli
Sequence	MKGWFLVLVIAIVGEVIATSALKSSEGFTKLAPSAVVIIGYGIIFYFLSLVLKSIPVG VAY AVWSGLGVVITAIWLLHGGKLDWGFVGMGLIVSGVVVLNLLSKASAH
Source	in vitro E.coli expression system
Target Names	qacE
Protein Names	Recommended name: Quaternary ammonium compound-resistance protein qacE Alternative name(s): Quaternary ammonium determinant E
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.