



Recombinant Human Mitochondrial inner membrane protein OXA1L (OXA1L)

Product Code	CSB-CF613583HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q15070
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAMGLMCGRRELLRLLQSGRRVHSVAGPSQWLKPLTTRLLFPVAPCCCRP HYLFLAASG PRSLSTSAISFAEVQVQAPPVVAATPSPTAVPEVASGETADVVTAAEQSFAE
Source	in vitro E.coli expression system
Target Names	OXA1L
Protein Names	Recommended name: Mitochondrial inner membrane protein OXA1L Alternative name(s): Hsa OXA1Hs Oxidase assembly 1-like protein Short name= OXA1-like protein
Expression Region	1-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.