



Recombinant Halobacterium salinarum Protein CrcB homolog 1 (crcB1)

Product Code	CSB-CF888117HTM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9HNW2
Product Type	Transmembrane Protein
Immunogen Species	Halobacterium salinarum (strain ATCC 700922 / JCM 11081 / NRC-1) (Halobacterium halobium)
Sequence	MTGAVAPPAVLVAAGGALGAVLRWRVVAATPTTEYPAGTLVVNVVGSFVLAA LTFAAADA DTMLLFGTGACGAFTTFASFSDVVALVDADRPVAAAGHALGNLLGAGLAVAL AWLLVA
Source	in vitro E.coli expression system
Target Names	crcB1
Protein Names	Recommended name: Protein CrcB homolog 1
Expression Region	1-119
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.