



Recombinant Escherichia coli O139:H28 Spermidine export protein MdtJ (mdtJ)

Product Code	CSB-CF421943EJD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZM59
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Sequence	MYIYWILLGLAIATEITGTLSMKWASVSEGNNGGFILMLVMISLSYIFLSFAVKKIAL GVA YALWEGIGILFITLFSVLLFDESLSLMKIAGLTTLVAGIVLIKSGTRKARKPELEVN HGA V
Source	in vitro E.coli expression system
Target Names	mdtJ
Protein Names	Recommended name: Spermidine export protein MdtJ
Expression Region	1-121
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.