



Recombinant Bacillus cereus Holin-like protein CidA (cidA)

Product Code	CSB-CF508381BQK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C1ENB1
Product Type	Transmembrane Protein
Immunogen Species	Bacillus cereus (strain 03BB102)
Sequence	MKWWKLSGQILLLFCFAWTGEWIAKQAHLPPVPGSIIGIFLLLISLKFNLVKKEWI QDGAD FLLKELILFFIPSAVAVIRYKDTLSQYGIDLILIMISTLCVTLVTGLLTELKLRKGS V Q
Source	in vitro E.coli expression system
Target Names	cidA
Protein Names	Recommended name: Holin-like protein CidA
Expression Region	1-121
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.