



Recombinant Rhodopseudomonas palustris Protein CrcB homolog (crcB)

Product Code	CSB-CF352780RLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61394
Product Type	Transmembrane Protein
Immunogen Species	Rhodopseudomonas palustris (strain ATCC BAA-98 / CGA009)
Sequence	MAYLLVFGGGLGAMFRHFINTLSGROLLGTAFPYHTFFINVTGSIVMGLIAGYLA FKGGS SQHFRLFLMTGILGGYTTFSAFSLDAALLYERGAVGLAVVYVLGSSVLAIALGLF GGMALI RAMT
Source	in vitro E.coli expression system
Target Names	crcB
Protein Names	Recommended name: Protein CrcB homolog
Expression Region	1-124
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.