



Recombinant *Desulfovibrio vulgaris* Protein CrcB homolog (crcB)

Product Code	CSB-CF491691DJR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B8DJU6
Product Type	Transmembrane Protein
Immunogen Species	<i>Desulfovibrio vulgaris</i> (strain Miyazaki F / DSM 19637)
Sequence	MQKIVLLGLAGALGSLARYGLAGLVQRAAPGSFPLGTFIVNVLGCLAFGFVWG VCENRIS LHPDLRVVLLTGFMGAFTTFSTFTFESLGLMETGQWLAFALYAGGQLLLGLALL WLGLGT GRLV
Source	in vitro E.coli expression system
Target Names	crcB
Protein Names	Recommended name: Protein CrcB homolog
Expression Region	1-124
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.