



Recombinant Mouse Protein tweety homolog 2 (Ttyh2)

Product Code	CSB-CF668837MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3TH73
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	TCAGPRAWKYFINRDRDYDDIDDDDPFNPQARRIAAHNPTRGQLHSFCSYSS GLGSQCSL QPPSQTISNAPVSEYMNQAILFGGNPRYENVPLIGRGSPPPTYSPSMRPTYMS VADEHLR HYEFPS
Source	in vitro E.coli expression system
Target Names	Ttyh2
Protein Names	Recommended name: Protein tweety homolog 2 Short name= mTTY2
Expression Region	407-532
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.