



Recombinant Chromobacterium violaceum Protein CrcB homolog (crcB)

Product Code	CSB-CF755858CKA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7NWP1
Product Type	Transmembrane Protein
Immunogen Species	Chromobacterium violaceum (strain ATCC 12472 / DSM 30191 / JCM 1249 / NBRC 12614 / NCIMB 9131 / NCTC 9757)
Sequence	MWKSILAI AIGAAAGALLRWFLGLRLNSLFPSPGTLAANLVGGYIIGVAVALF ADMPM LSPLWRLLIITGFCGGLTTFSTFSAEVVDLLRQGLQPAFAAIAVHVSQSLLMT MAGIAS WQWLRR
Source	in vitro E.coli expression system
Target Names	crcB
Protein Names	Recommended name: Protein CrcB homolog
Expression Region	1-126
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.