



Recombinant Escherichia coli O17:K52:H18 Protein CrcB homolog (crcB)

Product Code	CSB-CF483687EOG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7N9N0
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Sequence	MLQLLLAVFIGGGTGSVARWLLSMRFNPLHQAIPLGTLAANLIGAFIIGMGFAW FSRMTN IDPVWKVLITTFGCGGLTTFSTFSAEVVFLQEGRLGWALLNVFVNLLGSFAMT ALAFWL FSASTAH
Source	in vitro E.coli expression system
Target Names	crcB
Protein Names	Recommended name: Protein CrcB homolog
Expression Region	1-127
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.