



Recombinant Escherichia coli Protein CrcB (crcB)

Product Code	CSB-CF336246ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P37002
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MLQLLLAVFIGGGTGSVARWLLSMRFNPLHQAIPLGTLTANLIGAFIIGIGFAWF SRMTN IDPVWKVLITTGFCGGLTTFSTFSAEVVFLQEGRFGWALLNVFVNULLGSFAMT ALAFWL FSASTAH
Source	in vitro E.coli expression system
Target Names	crcB
Protein Names	Recommended name: Protein CrcB
Expression Region	1-127
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.