



Recombinant Escherichia coli Probable 4-amino-4-deoxy-L-arabinose-phosphoundecaprenol flippase subunit ArnF (arnF)

| | |
|--------------------------|---|
| Product Code | CSB-CF478194ENW |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B6I7K2 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Escherichia coli (strain SE11) |
| Sequence | MGLIWGLFSVIIASVAQLSLGFAASHLPPMTHLWDFIAALLAFGLDARILLGLLGYLLS VFCWYKTLHKLALSKAYALLSMSYVLVWIASMVLPGWEGTFSLKALLGVACIMSGLMLIF LPTTKQRY |
| Source | in vitro E.coli expression system |
| Target Names | arnF |
| Protein Names | Recommended name: Probable 4-amino-4-deoxy-L-arabinose-phosphoundecaprenol flippase subunit ArnF Short name= L-Ara4N-phosphoundecaprenol flippase subunit ArnF Alternative name(s): Undecaprenyl phosphate-aminoarabinose flippase subunit ArnF |
| Expression Region | 1-128 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |