



Recombinant Escherichia coli O17:K52:H18 Probable 4-amino-4-deoxy-L-arabinose- phosphoundecaprenol flippase subunit ArnF (arnF)

Product Code	CSB-CF482802EOG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7N5M4
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Sequence	MGLMWGLFSVIIASAAQLSLGFAASHLPPMTHLWDFIAALLAFGLDARILLGLL GYLLS VFCWYKTLHKLALSKAYALLSMSYVLVWIASMILPGWEGTFSLKALLGVACIMS GLMLIF LPTTKQRY
Source	in vitro E.coli expression system
Target Names	arnF
Protein Names	Recommended name: Probable 4-amino-4-deoxy-L-arabinose- phosphoundecaprenol flippase subunit ArnF Short name= L-Ara4N- phosphoundecaprenol flippase subunit ArnF Alternative name(s): Undecaprenyl phosphate-aminoarabinose flippase subunit ArnF
Expression Region	1-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.